

ResIWave 77

The professional **Resonance Initiated Wave** generator

- flexible application
- small and easy to use
- Med CE Certified in Europe
- Large range of accessories

Modern SMD-technology
with 4layer- PCB



Innovation

... in frequency therapy

according to Dr. Hulda Clark, Dr. Royal Rife and Dr. Robert Beck.

Chip-cards

with the latest chip-card technology:

- Freely programmable and from pre-programmed chip-cards
- Chip-cards according to Dr. Clark, Dr. Rife and Dr. Robert Beck
- Includes a number of built-in ready to run programs
- Resonance Initiated Frequency Effect (R.I.F.E.) class unit



Package contents:

- 1 *ResiWave* generator 77
- 1 Database chip-card with top 100 frequencies according to Clark, Rife, Beck and Baklayan
- 2 Red/black tear-proof long life cables with button connectors
- 2 Stainless steel hand electrodes (cylinder form)
- 2 Yellow wristbands*, non metallic
- 1 Battery*, 9 Volt E-Block
- 1 Portable case
- 1 User manual

* Consumable articles - no warranty

Configuration:

- Chip-card drive for upgrades/applications
- 3 x 7 min, with 20 minute break program
- Continuous 1 hour mode
- 1 KHz for zappicator use, plus additional programs for daily use based on research from Dr. Robert Beck and Alan E. Baklayan ND

30 programmable frequencies:

- Simple search by name from dBase chip-cards
- Individual setting of frequency: 1 Hz to 1 MHz
- Up to 30 frequencies can be freely programmed in internal memory and run automatically.
- Programs can be easily given a descriptive name and stored on reusable master chip-cards that can be easily run on the *ResiWave 77* and the *ResiWave 21* (playback only)
- Allows therapist to program an individual therapy for use on the *ResiWave 21* which the patient can use at home.
- The master chip-cards are the same size as SIM cards (like in cell phones) and can be reprogrammed an unlimited number of times.

Technical data:

- Frequency-range: 1 Hz to 1 MHz (.01 Hz steps)
- Frequency accuracy: 100 ppm (0.01%)
- Wave form: Square
- Offset: positive (Dr. Clark's standard), positive with 0.25 V remaining voltage no offset (plus/minus)
- Output voltage 1-15 volts in 0.2 volt steps, 1W
- Intensity can be adjusted during the program
- Frequency steps can be easily skipped/repeated
- Automatic contact monitoring (override possible)
- Softstart: Amplitude changes slowly to prevent an unpleasant tingling sensation on frequency change
- Display: 2 lines, 12 characters (letters, numbers and special characters) easy user guidance and parameter display. Remaining time during program
- 19 Key membrane keypad (for entering numbers and characters as well as menu navigation)
- Several operating languages available
- Automatic self-test after the unit is switched on
- Power supply: 9 V E-block battery (8.4 V NiMH rechargeable battery, 250 mA/h)
- Firmware updates possible via chip-card

Accessories: (not included with product)

Database chip-cards:

- Chip-card with approx. 350 frequencies according to Dr. Clark in groups and sorted alphabetically
- Chip-card with approx. 5000 frequencies according to Dr. Rife sorted alphabetically
- Chip-card with frequencies according to A.E. Baklayan ND, etc.

Chip-cards:

- Range of pre-programmed chip-cards for use on *ResiWave 21* or *77* can be used (new chip-cards are constantly being released, see website for details)
- Master chip-card: Stores up to 30 frequencies for use on *ResiWave 21* and *77*. Can be reprogrammed any number of times on the *ResiWave 77* only

Other Accessories:

Zappicator, Plate Zapping, Gold Plates, Brain-tuner upgrade, Colloidal Silver Upgrade, Carbon electrodes, Stainless steel electrodes


More details:

You can learn more about the *ResiWave 77* (for practitioners) and the *ResiWave 21* (for users) on our website:



<http://www.resiwave.com>

Distributor:



Disclaimer:

The following information is for therapists only! Reports on the effectiveness and beneficial results from all *ResiWave* products are based on user accounts, rather than acceptance by conventional medicine or science in general. These instruments are sold as electronic test instruments. No claims regarding medical diagnoses, treatments or cures are stated or implied. The BW21 and BW77 units are approved as medical devices for use in Europe (med CE, class IIa, rule 9). They have not been evaluated by the FDA.